



# High Content Imaging



## SlideExpress 2

Microscope Automated  
Slide Imaging System



# Key Features

## The Perfect Platform for High Content Slide Imaging

An Ni-E motorized upright microscope platform with several choices for high performance optics and a wide range of available detectors means the highest quality results.

Using the Marzhauser Slide Express 2 automated slide loader, stability and repeatability is maintained by a solid platform with precision connection to the microscope body, resulting in flawless loading and unloading of samples.

With a unified software package across all microscope platforms, data can be easily shared and visualized as well as learning time reduced.



# STABLE AND RELIABLE



Slides are inserted into the cassette using two-slide holders, which insure the slides are flat and oriented properly for imaging.

Cassettes are loaded into a motorized enclosure with space for up to 3 cassettes, or a total of 120 slides.

To place slides on the imaging area, the two slide holders are extracted from the cassettes by magnetic attachment. Because the stage is integrated and attached to the loading mechanism, a stable and reliable transfer of slides results.



# NIS-Elements Software






NIS-Elements reads and catalogs the cassettes, and automates the delivery of specimens to the microscope. An indicator shows which cassette location is currently loaded onto the microscope stage.

Slide scanning can be in a defined area and/or pattern, or specimen boundaries can be detected and the area within used as a tiling, fixed, or random pattern.


These areas can be defined at the time each slide is loaded, or a predefined area used for all slides. Users

can select specific slides to acquire, or choose to scan all slides in all cassettes in one run.


Imaging can be performed in brightfield and/or epifluorescence modalities, with results tiled or presented as image stacks.



# Acquisition




Slides can be quickly scanned at low magnification to locate items of interest, including several separate items per slide...



...then rescanned at high resolution and high magnification, so that only the areas of the slide with useful data are acquired and tiled.

High resolution images make up the tiled scan. Not acquiring blank areas results in much faster scanning and less storage space required per slide.



# Image Analysis





NIS-Elements' **General Analysis** engine employs a powerful processing and analysis toolbox for users.

Dedicated turnkey assays can be applied to image data as well as user-customized assays for specific applications.



All routines can be executed during acquisition, for example: to modify the course of an experiment, or post-acquisition on the dedicated server.

Assay results can be collated and statistical information displayed, or be exported.



Some examples are shown here.



There are many tools which can be customized to detect and automatically segment and classify objects. In this example, pituitary gland sections were imaged over a large field of view, tiled, and colloids were auto-detected and measured.




This example tool computes areas occupied by different fluorescence markers in tiled brain images.



This tissue microarray punch was segmented using an IHC detection tool which automatically finds hematoxylin and diaminobenzidine (DAB) colors.


# Specifications



## Microscope Specifications

<b>Imaging Modalities</b>	Transmitted Light Brightfield Widefield Epifluorescence Point Scanning confocal (option with C2)
<b>Illumination</b>	8 channel Epifluorescence LED LED transmitted Light Laser illumination for confocal sources (option)
<b>Objectives</b>	All air objectives lenses can be used for high content
<b>Imaging Methods</b>	Multidimensional XY, Z, wavelength, multi-stage position including tiling
<b>Autofocus</b>	Software contrast-based autofocus
<b>Hardware Triggering</b>	Supports direct hardware triggering of light sources and shutters
<b>Antivibration Table</b>	Included 30x48" x 6" table with full system
<b>Dimensions (without detector*)</b>	32cm x 60cm x 510cm typical (height increases 70mm with additional EPI port)
<b>Power Consumption</b>	

*\*see detector brochures for detector specification details*



## Slide Loader Specifications

<b>Device</b>	Marzhauser Wetzlar SlideExpress2
<b>Application</b>	Slide loading
<b>Capacity</b>	3 Cassettes each holding up to 40 1x3 inch glass slides: total of 120 slides
<b>Robotics</b>	Internal elevator robot and picker robot
<b>Communication</b>	USB through Tango controller
<b>Power Requirements</b>	100V / 50Hz 115V / 60Hz / 240V / 50Hz ; switchable
<b>Stage Specifications</b>	Accuracy +/- 1um; repeatability <1um bidirectional; resolution 0.05um step; max speed 100mm/s
<b>Power Consumption</b>	<100W
<b>Dimensions</b>	108 x 75 x 108.5cm

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer.  
April 2020 ©2020 Nikon Instruments B.V.

 <b>WARNING</b>	TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.
--	--

Company names and product names appearing in this brochure are their registered trademarks or trademarks.

<b>WARNING-LASER RADIATION AVOID EXPOSURE TO BEAM CLASS 3B LASER PRODUCT</b> Total Power 500mW MAX. CW 400~700nm IEC/EN60825-1 : 2007, 2014
Complies with FDA performance standards for laser products except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007

<b>DANGER-VISIBLE AND INVISIBLE LASER RADIATION AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION CLASS 4 LASER PRODUCT</b> Total Power 1500mW MAX. CW 370~790nm IEC/EN60825-1 : 2007, 2014
Complies with FDA performance standards for laser products except for deviations pursuant to Laser Notice No.50 dated June 24, 2007.



**NIKON INSTRUMENTS EUROPE BV.**  
Tripolis 100, Burgerweeshuispad 101, 1076 ER Amsterdam, Netherlands  
Enquiries Tel: 0031-(0)-20-7099-000

[www.microscope.healthcare.nikon.com](http://www.microscope.healthcare.nikon.com)